

From Page No. 15

Want to try a serum free (PS04) run
on CHO cells to check for Fus expression
levels

Have 3 sets of plates currently plus
6 well dishes

Replaced media on 6 well dishes
w/ PS04 → Inc

FOR CHO → needs insulin (bovine)
transferrin
Lipid
trace elements
HEPES

Inc 37°C →

To Page No. 17

Witnessed & Understood by me,

Date

Invented by

Date THURS

Recorded by

6/17/93

Project No. 1713
Book No. 18002

17

TITLE _____


From Page No. 16

Exhibit H, pg. 2 of 20

Froze down 1 set of FUS picks in
serum + 7% DMSO -70°C.

Split 10cm plates each 1:10.

Cont P504 on 6 well dishes



To Page No. 18

Witnessed & Understood by me, _____

Date _____

Invented by _____

Date FRI

Recorded by _____

W. M. Bacon

6/18/93

Project No. 1713
Book No. 18002TITLE From Page No. 17

Harvested P504 media's
Filtered each Through syringe 0.45µm
Stored -20°C o/n

Split all 36 Fus picks.

Witnessed & Understood by me,

Date

Invented by

Date MON

Recorded by

6/21/93

Project No. 1713

Book No. 18002

19


TITLE _____

From Page No. 12

Dot Blotted 200µl each P504 media
Decorated blot w/ 2 Human Fc
ELL Detected

#'s 9, 11 + 22 appear to have highest
expression → discarded all cultures
except 9, 11 + 22.

Ready to scale-up for Protein production



Witnessed & Understood by me.

Date

Invented by

Date

Recorded by

6/22/93

Will Baron

9

1713

18002


TITLE

From Page No 12

Dot Blotted 200 µl each P304 media
Decorated blot w/ 2 Human Fc
ELL Detected

#'s 9, 11 + 22 appear to have highest
expression → discarded all cultures
except 9, 11 + 22.

Ready to scale-up for Protein production



Witnessed & Understood by me.	Date	Invented by	Date
		Recorded by	

WILL BAYON 6/22/93

Project No. 1713
Book No. 18002 TITLE _____

20

From Page No. 19

Split all FNS samples 1-10 in 10cm dishes
(9, 11 + 22)

& 1-5 on 15cm dishes (~~#5~~ #9 only)

To Page No. 21

Witnessed & Understood by me,

Date

Invented by

Date FR1

Recorded by

WMP
Baton
6/25/93

Project No. 1713
Book No. 18002

TITLE _____

From Page No. 20

Split Fns 11 & 22 each 1:10 on 10cm dishes.
Split Fns 9 15cm plate to 5x 15cm plates
Getting ready for 1st P504 run.

To Page No. 21

Witnessed & Understood by me, _____

Date _____

Invented by _____

Recorded by _____

Date MON6/28/93

From Page No. 21

Split each Fus 9 15cm plate to 8x15cm plates
(= 5x8 = 40 total)

Inc o/n 37°C 5% CO₂

To Page No. 23

Witnessed & Understood by me,

Date

Invented by

Date WED

Recorded by

6/30/93

111

Projec. 5. 1713
Book No. 18002

Exhibit H, pg. 9 of 20

23

TITLE _____

From Page No. 22

Changed media on 40 x 15 cm FUS 9 plates
to PSO4 → Inc 37°C

Will harvest in ~ 4-6 days (when cells
lift off)

Started FUS 9 in 250ml spinner
flask (with 1 entire 16cm dish
worth of cells).

To Page No. 24

Witnessed & Understood by me,

Date

Invented by

Date 5/4/93

Recorded by

W. M. B. B. B.
7/3/93

Project No. 713
Book No. 18002 TITLE _____

-24

From Page No. 23

Split Fus 9, 11 + 22 all 1:10

Checked Fus 9 spinner → split 1:10

FS04 cont until 7/7

To Page No. 25

Witnessed & Understood by me,

Date

Invented by

Recorded by

Date TUES

7/6/93

TITLE _____


From Page No. 29

Harvested PS04 media

Filtered through 0.45 μ m filter

Added protease inhibitors (PMSF, aprotinin, leupeptin
& pepstatin)

Stored 4°C o/n.



To Page No. 26

Witnessed & Understood by me,

Date

Invented by

Date WED

Recorded by

7/7/93

From Page No. 25Protein A column procedure:

Column is washed w/ 0.1M Na Citrate pH 6.0
Sample is made up to 0.1M Na Citrate pH 6.0
filtered & then loaded @ ~4ml/min
after loading column is washed w/ citrate
→ when baseline level is reached
then sample is eluted w/ MgCl₂
elution buffer (see lab protocols)
Sample is desalted on a PD-10 column
& concentrated w/ 545 unit.

Ran entire 1st P504 batch, washed,
eluted, desalted & stored 4°C o/n
(Volume is ~3.5ml in PBS)

To Page No. 27

Witnessed & Understood by me,

Date

Invented by

Recorded by

Date

THURS

7/8/93

TITLE _____


Exhibit H, pg. 13 of 20

From Page No. 26

Ran samples of Pre Prot A, Post Prot A & of recovered material on 10% SDS-PAGE in duplicate for analysis

MW	MW	50µl REC	25µl REC	10µl REC	5µl POST	50µl PRE	MW	MW	MW	50µl REC	25µl REC	10µl REC	50µl POST	50µl PRE	MW
16	15	14	13	12	11	10	9	8	7	6	5	4	3	2	1

Cut gel into 2 halves → WESTERN transferred one half & Ponceauid → other half coomassie stained



Did BCA (on 100µl of 3.5ml) for protein concentration determination

A₅₆₂ recovered FUS 9 = 0.192
(see standard curve p. 28)

~~0.192 x 29.27 = 5.62094~~

$$y = 29.27(0.192) - 0.60091$$

$$y = 5.02 \mu\text{g in } 100\mu\text{l}$$

$$= 50\text{ng}/\mu\text{l}$$

$$= 175\mu\text{g total yield!}$$

To Page No. 28


Invented by

Date FRI

Recorded by

7/9/93

27



I think That I will get much better
yields from spinner culture.

Split FUS 9, 11 & 22 plates all 1:10

Split FUS 9 spinner 1:10

Witnessed & Understood by me.

Date


Invented by

Date

Recorded

W. M. Bacon

F21
7/9/03



TITLE _____

Exhibit H, pg. 16 of 20

From Page No. 28

Decorated FUS 9 western w/ α Human Fc
ECL detected - Many species observed
including some probable degradation products
Photographed coomassie stained portion of gel

↓ here major material
is seen at $\geq 200K$

Amly
11
50 μ l
25 μ l
10 μ l
POST
PRE
11
Amly

FUS 9
WESTERN

Split all cells & spinner

To Page No. 3

Witnessed & Understood by me,

Date

Invented by


Recorded by

Date MON

7/12/93

From Page No. 33

Stopped o/n gel → cut into 2 halves

Coomassie stained $\frac{1}{2}$
WESTERN transferred other half → Ponceau stainedS.T
FA/RT/F

Stored blot RT

Destained Coomassie o/n

Started PS04 run #4 as before

To Page No. 35

Witnessed & Understood by me,

Date

Invented by

Recorded by


WIM BAYON

Date FRI

7/23/93

35

red



272193
15101P


2439-76

TITLE

Exhibit H, pg. 19 of 20

From Page No. 34

Decorated Fusion Western w/ α human Fc \rightarrow ECL detected
Photographed Coomassie stain



FUS
IgG
ECL

To Page No. 36

Witnessed & Understood by me,

Date

Invented by

Recorded by

Date MON

7/26/93

From Page No. 35

Harvested 1st P504 run
Ran prot A column
Eluted w/ $MgCl_2$
Desalted, stored 40c

Split all FUS plates + spinners (NOT P504's)

To Page No. 37

Witnessed & Understood by me,

Date

Invented by

Recorded by

Date MON7/26/93