

From Page No. 15

Want to try a serum free (PS04) run  
on CHO cells to check for Fus expression  
levels

Have 3 sets of plates currently plus  
6 well dishes

Replaced media on 6 well dishes  
w/ PS04 → Inc

FOR CHO → needs insulin (bovine)  
transferrin  
Lipid  
trace elements  
HEPES

Inc 37°C →

To Page No. 17

Witnessed &amp; Understood by me,

Date

Invented by

Date THURS

Recorded by

6/17/93

Project No. 1713  
Book No. 18002

17

TITLE \_\_\_\_\_


From Page No. 16

Exhibit H, pg. 2 of 20

Froze down 1 set of FUS picks in  
serum + 7% DMSO -70°C.

Split 10cm plates each 1:10.

Cont P504 on 6 well dishes



To Page No. 18

Witnessed & Understood by me, \_\_\_\_\_

Date \_\_\_\_\_

Invented by \_\_\_\_\_

Date FRI

Recorded by \_\_\_\_\_

W. M. Bacon

6/18/93

From Page No. 17

Harvested PSCY media's  
Filtered each Through syringe 0.45µm  
Stored -20°C o/N

Split all 36 Fus picks.

Witnessed &amp; Understood by me,

Date

Invented by

Date MON

Recorded by

6/21/93

Project No. 1713

Book No. 18002

19


TITLE \_\_\_\_\_

From Page No. 12

Dot Blotted 200µl each P504 media  
Decorated blot w/ 2 Human Fc  
ELL Detected

#'s 9, 11 + 22 appear to have highest  
expression → discarded all cultures  
except 9, 11 + 22.

Ready to scale-up for Protein production



Witnessed & Understood by me.

Date

Invented by

Date

Recorded by

6/22/93

*Will Baron*

9

1713

18002


FILE

From Page No 12

Dot Blotted 200 µl each P304 media  
Decorated blot w/ 2 Human Fc  
ELL Detected

#'s 9, 11 + 22 appear to have highest  
expression → discarded all cultures  
except 9, 11 + 22.

Ready to scale-up for Protein production



Witnessed & Understood by me.

Date

Invented by

Recorded by

Will Brown

Date TUES

6/22/93

Project No. 1713  
Book No. 18002 TITLE \_\_\_\_\_

20

From Page No. 19

Split all FNS samples 1-10 in 10cm dishes  
(9, 11 + 22)

& 1-5 on 15cm dishes (~~#5~~ #9 only)

To Page No. 21

Witnessed &amp; Understood by me,

Date

Invented by

Date FR1

Recorded by

WMP  
Bacon  
6/25/93

Project No. 1713  
Book No. 18002

TITLE \_\_\_\_\_

From Page No. 20

Split Fns 11 & 22 each 1:10 on 10cm dishes.  
Split Fns 9 15cm plate to 5x 15cm plates  
Getting ready for 1st P504 run.

To Page No. 21

Witnessed &amp; Understood by me, \_\_\_\_\_

Date \_\_\_\_\_

Invented by \_\_\_\_\_

Recorded by \_\_\_\_\_

Date MON6/28/93

From Page No. 21

Split each Fus 9 15cm plate to 8x15cm plates  
(=  $5 \times 8 = 40$  total)

Inc o/n 37°C 5% CO<sub>2</sub>

To Page No. 23

Witnessed &amp; Understood by me,

Date

Invented by

Date WED

Recorded by

6/30/93



111

Projec.

J. 1713

Exhibit H, pg. 9 of 20

Book No. 18002

23

TITLE \_\_\_\_\_

From Page No. 22

Changed media on 40 x 15 cm FUS 9 plates  
to PS04 → Inc 37°C

Will harvest in ~ 4-6 days (when cells  
lift off)

Started FUS 9 in 250ml spinner  
flask (with 1 entire 16cm dish  
worth of cells).

To Page No. 24

Witnessed &amp; Understood by me,

Date

Invented by

Date SAT

Recorded by

7/3/93

Project No. 713  
Book No. 18002 TITLE \_\_\_\_\_

-24

From Page No. 23

Split Fus 9, 11 + 22 all 1:10

Checked Fus 9 spinner → split 1:10

FS04 cont until 7/7

To Page No. 25

Witnessed &amp; Understood by me,

Date

Invented by

Recorded by

Date TUES

7/6/93

TITLE \_\_\_\_\_


From Page No. 29

Harvested PS04 media

Filtered through 0.45  $\mu$ m filter

Added protease inhibitors (PMSF, aprotinin, leupeptin  
& pepstatin)

Stored 4°C o/n.



To Page No. 26

Witnessed & Understood by me,

Date

Invented by

Date WED

Recorded by

W. P. Bacon  
7/7/93

From Page No. 25Protein A column procedure:

Column is washed w/ 0.1M Na Citrate pH 6.0  
Sample is made up to 0.1M Na Citrate pH 6.0  
filtered & then loaded @ ~4ml/min  
after loading column is washed w/ citrate  
→ when baseline level is reached  
then sample is eluted w/ MgCl<sub>2</sub>  
elution buffer (see lab protocols)  
Sample is desalted on a PD-10 column  
& concentrated w/ 545 unit.

Ran entire 1<sup>st</sup> P504 batch, washed,  
eluted, desalted & stored 4°C o/n  
(Volume is ~3.5ml in PBS)

To Page No. 27

Witnessed &amp; Understood by me,

Date

Invented by

Recorded by

Date

THURS

7/8/93

TITLE \_\_\_\_\_


Exhibit H, pg. 13 of 20

From Page No. 26

Ran samples of Pre Prot A, Post Prot A & of recovered material on 10% SDS-PAGE in duplicate for analysis

MW	MW	50µl REC	25µl REC	10µl REC	5µl POST	50µl PRE	MW	MW	MW	50µl REC	25µl REC	10µl REC	50µl POST	50µl PRE	MW
16	15	14	13	12	11	10	9	8	7	6	5	4	3	2	1

Cut gel into 2 halves → WESTERN transferred one half & Ponceauid → other half coomassie stained



Did BCA (on 100µl of 3.5ml) for protein concentration determination

A<sub>562</sub> recovered FUS 9 = 0.192  
(see standard curve p. 28)

~~0.192 x 29.27 = 5.6091~~

$$y = 29.27(0.192) - 0.60091$$

$$y = 5.02 \mu\text{g in } 100\mu\text{l}$$


$$= 50\text{ng}/\mu\text{l}$$

$$= 175\mu\text{g total yield!}$$

To Page No. 28

Invented by	Date
W. J. Brown	FRI 7/9/93
Recorded by	

27



I think That I will get much better  
yields from spinner culture.

Split FUS 9, 11 & 22 plates all 1:10

Split FUS 9 spinner 1:10

Witnessed &amp; Understood by me.


Date

Invented by

Date F21

Recorded

7/9/03



TITLE \_\_\_\_\_

Exhibit H, pg. 16 of 20

From Page No. 28

Decorated FUS 9 western w/  $\alpha$  Human Fc  
ECL detected - Many species observed  
including some probable degradation products  
Photographed coomassie stained portion of gel

Amly  
11  
50kD  
25kD  
10kD  
POST  
PRE  
11

here major material  
is seen at  $\geq 200K$

FUS 9  
WESTERN

Split all cells & spinner

To Page No. 3

Witnessed & Understood by me,

Date

Invented by

Date MON

Recorded by


7/12/93

W. M. Barton



From Page No. 33

Stopped o/n gel → cut into 2 halves

Coomassie stained  $\frac{1}{2}$   
WESTERN transferred other half → Ponceau stainedS.T  
FAL/F

Stored blot RT

Destained Coomassie o/n

Started PS04 run #4 as before

To Page No. 35

Witnessed &amp; Understood by me, \_\_\_\_\_

Date \_\_\_\_\_

Invented by \_\_\_\_\_

Recorded by \_\_\_\_\_


WIM BAYON

Date FRI

7/23/93

35

ted



2/2/93  
15:00P


2/2/93

TITLE

Exhibit H, pg. 19 of 20

From Page No. 34

Decorated Fusion Western w/  $\alpha$  human Fc  $\rightarrow$  ECL detected  
Photographed Coomassie stain



FUS  
IgG  
ECL

To Page No. 36

Witnessed & Understood by me,

Date

Invented by

Recorded by

Date MON

7/26/93

W. H. Bacon

Project No. 1713

Book No. 18002 TITLE \_\_\_\_\_

36

From Page No. 35

Harvested 1st P504 run  
 Run prot A column  
 Eluted w/  $MgCl_2$   
 Desalted, stored 40c

Split all PUS plates + spinners (NOT P504's)

To Page No. 37

Witnessed & Understood by me,

Date

Invented by

Date MON

Recorded by

WMB

7/26/93